

Omnion 3.0 Software Recommended Compatibiity Matrix

	Windows 98	Windows NT	Windows 2000	Windows XP	Windows Vista	Windows 7
QuikChem 8000						
3.0.218	X	X	X	X		
3.0.219	X	X	X	X		
3.0.220	X	X	X	X		
3.0.222				X		
3.0.223				X		
3.0.224*				X		
QuikChem 8500						
3.0.219	X	X	X	X		
3.0.220	X	X	X	X		
3.0.222				X		
3.0.223				X		
3.0.224*				X		
QuikChem 8500 Series 2						
3.0.222				X		
3.0.223				X	X	
3.0.224*				X	X	X

*release July 2010